

# XDOP Device Class Description

Document Version 1.0

Includes Descriptions for following devices:

- urn:schemas-coscom-org:device:MCU5Ladder:1

creation date: Thursday, September 24, 2009

## Table of Contents

|   |    |
|---|----|
| 1 Device class MCU5Ladder.....                    | 5  |
| 1.1 Service class SunTechTangoStressBPBridge..... | 7  |
| 1.1.1 Actions.....                                | 8  |
| 1.1.1.1 BeginBPStudy.....                         | 8  |
| 1.1.1.2 ClearBPReadings.....                      | 8  |
| 1.1.1.3 EndBPStudy.....                           | 8  |
| 1.1.1.4 RecallBPValues.....                       | 8  |
| 1.1.1.5 StartBPSample.....                        | 8  |
| 1.1.1.6 StopBPSample.....                         | 9  |
| 1.2 Service class TreadmillProfileEditor.....     | 10 |
| 1.2.1 Actions.....                                | 12 |
| 1.2.1.1 GetProfileHeader.....                     | 12 |
| 1.2.1.2 ReadProfileStep.....                      | 13 |
| 1.2.1.3 WriteUserProfileStep.....                 | 14 |
| 1.3 Service class TreadmillGeneralConfig.....     | 16 |
| 1.3.1 Actions.....                                | 18 |
| 1.3.1.1 GetDefaultType.....                       | 18 |
| 1.3.1.2 GetErrors.....                            | 18 |
| 1.3.1.3 GetFirmwareVersion.....                   | 19 |
| 1.3.1.4 GetLockMask.....                          | 19 |
| 1.3.1.5 GetOEMCode.....                           | 20 |
| 1.3.1.6 GetRTCTime.....                           | 20 |
| 1.3.1.7 GetSerialNumber.....                      | 21 |
| 1.3.1.8 GetType.....                              | 21 |
| 1.3.1.9 SetFailSafeTimeout.....                   | 21 |
| 1.3.1.10 SetProtocolParams.....                   | 22 |
| 1.3.1.11 SetRTCTime.....                          | 23 |
| 1.3.1.12 SetSerialNumber.....                     | 23 |
| 1.4 Service class TreadmillUserControl.....       | 25 |
| 1.4.1 Actions.....                                | 29 |
| 1.4.1.1 Cooldown.....                             | 29 |
| 1.4.1.2 Countdown.....                            | 29 |
| 1.4.1.3 DecGrade.....                             | 30 |
| 1.4.1.4 DecSpeed.....                             | 30 |
| 1.4.1.5 FreezeStepperTestPause.....               | 30 |
| 1.4.1.6 GetHRCStatus.....                         | 30 |
| 1.4.1.7 GetNextStepStatus.....                    | 31 |
| 1.4.1.8 GetPersonData.....                        | 32 |
| 1.4.1.9 GetProfileStatus.....                     | 32 |
| 1.4.1.10 GetProfListEntry.....                    | 33 |
| 1.4.1.11 GetProfListSize.....                     | 34 |
| 1.4.1.12 GetTestListEntry.....                    | 34 |
| 1.4.1.13 GetTestListSize.....                     | 34 |
| 1.4.1.14 GetTestStatus.....                       | 35 |
| 1.4.1.15 GetUKKTestResult.....                    | 35 |
| 1.4.1.16 HoldGrade.....                           | 36 |
| 1.4.1.17 HoldSpeed.....                           | 36 |
| 1.4.1.18 IncGrade.....                            | 37 |
| 1.4.1.19 IncSpeed.....                            | 37 |
| 1.4.1.20 Pause.....                               | 37 |
| 1.4.1.21 ProfileStepBack.....                     | 37 |
| 1.4.1.22 ProfileStepForward.....                  | 38 |
| 1.4.1.23 QuickStop.....                           | 38 |
| 1.4.1.24 SetPersonData.....                       | 38 |
| 1.4.1.25 SetProfileStepInterval.....              | 39 |

|   |    |
|---|----|
| 1.4.1.26 ShowNextStep.....                    | 39 |
| 1.4.1.27 StartHeartRateControl.....           | 40 |
| 1.4.1.28 StartManual.....                     | 40 |
| 1.4.1.29 StartNextStepperTestInterval.....    | 40 |
| 1.4.1.30 StartProfile.....                    | 41 |
| 1.4.1.31 StartTestConconi.....                | 42 |
| 1.4.1.32 StartTestCooper.....                 | 42 |
| 1.4.1.33 StartTestGardner.....                | 42 |
| 1.4.1.34 StartTestProfile.....                | 43 |
| 1.4.1.35 StartTestRamp.....                   | 43 |
| 1.4.1.36 StartTestStepper.....                | 43 |
| 1.4.1.37 StartTestUKK.....                    | 44 |
| 1.4.1.38 Stop.....                            | 45 |
| 1.4.1.39 TestCooperStepBack.....              | 45 |
| 1.4.1.40 TestCooperStepForward.....           | 45 |
| 1.5 Service class TreadmillOptions.....       | 46 |
| 1.5.1 Actions.....                            | 47 |
| 1.5.1.1 GetListEntry.....                     | 47 |
| 1.5.1.2 GetListSize.....                      | 47 |
| 1.5.1.3 GetOption.....                        | 48 |
| 1.5.1.4 GetVersion.....                       | 49 |
| 1.5.1.5 SetOption.....                        | 49 |
| 1.6 Service class TreadmillDirectControl..... | 51 |
| 1.6.1 Actions.....                            | 53 |
| 1.6.1.1 GetDriveAccelDecelRange.....          | 53 |
| 1.6.1.2 GetDriveSpeedRange.....               | 53 |
| 1.6.1.3 GetDriveStatus.....                   | 53 |
| 1.6.1.4 GetElevatorGradeRange.....            | 54 |
| 1.6.1.5 GetElevatorStatus.....                | 54 |
| 1.6.1.6 SetDriveSpeed.....                    | 55 |
| 1.6.1.7 SetElevatorGrade.....                 | 55 |
| 1.7 Service class TerminalLCD6x4Panel.....    | 57 |
| 1.7.1 Actions.....                            | 61 |
| 1.7.1.1 Beep.....                             | 61 |
| 1.7.1.2 GetDistanceUnit.....                  | 61 |
| 1.7.1.3 GetGradeUnit.....                     | 61 |
| 1.7.1.4 GetKeyLock.....                       | 61 |
| 1.7.1.5 GetSpeedUnit.....                     | 62 |
| 1.7.1.6 KeyDown.....                          | 62 |
| 1.7.1.7 KeyUp.....                            | 62 |
| 1.7.1.8 SetDistanceUnit.....                  | 63 |
| 1.7.1.9 SetGradeUnit.....                     | 63 |
| 1.7.1.10 SetKeyLock.....                      | 63 |
| 1.7.1.11 SetModeLED.....                      | 64 |
| 1.7.1.12 SetSpeedUnit.....                    | 64 |
| 1.8 Service class TreadmillMeasures.....      | 65 |
| 1.8.1 Actions.....                            | 67 |
| 1.8.1.1 ResetDistance.....                    | 67 |
| 1.8.1.2 ResetEnergy.....                      | 67 |
| 1.8.1.3 ResetTime.....                        | 67 |
| 1.8.1.4 SetDistance.....                      | 67 |
| 1.8.1.5 SetEventMask.....                     | 67 |
| 1.8.1.6 SetTime.....                          | 68 |
| 1.9 Service class MCU5FlashUpdate.....        | 69 |
| 1.9.1 Actions.....                            | 71 |
| 1.9.1.1 GetDefaultType.....                   | 71 |
| 1.9.1.2 GetFirmwareVersion.....               | 71 |
| 1.9.1.3 ReadBuffer.....                       | 71 |

|                                    |    |
|------------------------------------|----|
| 1.9.1.4 ShowUpdateScreen.....      | 72 |
| 1.9.1.5 SoftReset.....             | 72 |
| 1.9.1.6 StartUpdate.....           | 72 |
| 1.9.1.7 WriteBuffer.....           | 73 |
| 1.10 Service class DataLogger..... | 74 |
| 1.10.1 Actions.....                | 75 |
| 1.10.1.1 SetColumnDelimiter.....   | 75 |
| 1.10.1.2 SetRecordColumns.....     | 75 |
| 1.10.1.3 SetRecordEventPeriod..... | 75 |
| 1.10.1.4 StartLogging.....         | 76 |
| 1.10.1.5 StopLogging.....          | 76 |

## 1 Device class MCU5Ladder

Namespace: XDOP.schemas\_coscom\_org.device  
Assembly: MCU5Ladder.dll

### Syntax

```
public class MCU5Ladder : IDeviceObject
```

### Constructors

|              |  |
|--------------|--|
| MCU5Ladder() | Initializes a new instance of the class. |
|--------------|--|

### Device Properties

| Property name    | Value                                      |
|------------------|--|
| DeviceType       | urn:schemas-coscom-org:device:MCU5Ladder:1 |
| EPC              |  |
| FriendlyName     | Ladder ergometer                           |
| Manufacturer     | h/p/cosmos                                 |
| ManufacturerURL  | www.h-p-cosmos.com                         |
| ModelDescription |  |
| ModelName        |  |
| ModelNumber      |  |
| ModelURL         |  |
| SerialNumber     |  |

### Common Properties

| Property name       | Description  |
|---------------------|--|
| AcceptNewerVersions | Indicates if the device object accepts newer versions of the device type.                  |
| LazyLoading         | If set, the device loads its services not until they are really needed.                    |
| XDOPDevice          | Gets the XDOPLib.IDevice Interface of the device object if already connected               |
| Services            | Get's a read-only collection of the device services. For a list of services see 'Services' |
| ConnectionStatus    | Gets the actual XDOPLib.ConnectionStatus object for this device.                           |

### Device Services

Following list enumerates the service (type) objects of the 'Services' Property.

- [urn:schemas-coscom-org:service:SunTechTangoStressBPBridge:1](#)
- [urn:schemas-coscom-org:service:TreadmillProfileEditor:1](#)
- [urn:schemas-coscom-org:service:TreadmillGeneralConfig:1](#)
- [urn:schemas-coscom-org:service:TreadmillUserControl:1](#)
- [urn:schemas-coscom-org:service:TreadmillOptions:1](#)
- [urn:schemas-coscom-org:service:TreadmillDirectControl:1](#)
- [urn:schemas-coscom-org:service:TerminalLCD6x4Panel:1](#)
- [urn:schemas-coscom-org:service:TreadmillMeasures:1](#)
- [urn:schemas-coscom-org:service:MCU5FlashUpdate:1](#)
- [urn:schemas-coscom-org:service:DataLogger:1](#)

## Common Methods

| Name   | Description   |
|--|---|
| <code>void Connect(ITransportObject transportObject);</code> | Connects to a device.                                       |
| <code>void Disconnect();</code>                              | Closes a established the connection.                        |
| <code>ProtocolParameter GetProtocolParameter();</code>       | Gets the protocol parameter of the current connection.      |
| <code>ITransportObject GetTransportObject();</code>          | Gets the ITransportObject object of the current connection. |
| <code>void DisposeDevice();</code>                           | Disposes the current device object.                         |

## 1.1 Service class SunTechTangoStressBPBridge

This service manages a communication bridge to the SunTech Tango Blood Pressure device. If the Tango device is connected with one of the other serial ports, it can be controlled with the actions of this service.

See SunTech Tango manual for more details or contact [h/p/cosmos](mailto:h/p/cosmos).

Namespace: XDOP.schemas\_coscom\_org.service

Assembly: SunTechTangoStressBPBridge.dll

### Syntax

```
public class SunTechTangoStressBPBridge : IServiceObject
```

### Constructors

|  |  |
|--|--|
| SunTechTangoStressBPBridge(XDOPLib.IService service) | Initializes a new instance of the class. |
|--|--|

### Common Service Properties

| Property name | Value   |
|---------------|---|
| ServiceId     | urn:schemas-coscom-org:serviceId:TangoBP                    |
| MajorVersion  | 1   |
| MinorVersion  | 1   |
| ServiceType   | urn:schemas-coscom-org:service:SunTechTangoStressBPBridge:1 |

### Service Properties and Events

| Type          | Name    | Flags | Description   |
|---------------|---------|-------|---|
| System.String | Message | E G   | Message string from the Tango blood pressure device |

### Service Methods

| Name                            | Description                    |
|---------------------------------|--------------------------------|
| <a href="#">BeginBPStudy</a>    | Begin blood pressure study.    |
| <a href="#">ClearBPReadings</a> | Clear blood pressure readings. |
| <a href="#">EndBPStudy</a>      | End blood pressure study.      |
| <a href="#">RecallBPValues</a>  | Recall blood pressure values.  |
| <a href="#">StartBPSample</a>   | Start blood pressure sample.   |
| <a href="#">StopBPSample</a>    | Stop blood pressure sample.    |

### Service error codes

| Error code | Error description     |
|------------|-----------------------|
| 823        | No SunTech BP monitor |

### Details

On the following pages you will find a detailed description for every action and variable of this service.

### 1.1.1 Actions

#### 1.1.1.1 Action BeginBPStudy

Begin blood pressure study.

#### Syntax

```
public void BeginBPStudy()
```

#### Action Errors

|     |                       |
|-----|-----------------------|
| 823 | No SunTech BP monitor |
|-----|-----------------------|

#### 1.1.1.2 Action ClearBPReadings

Clear blood pressure readings.

#### Syntax

```
public void ClearBPReadings()
```

#### Action Errors

|     |                       |
|-----|-----------------------|
| 823 | No SunTech BP monitor |
|-----|-----------------------|

#### 1.1.1.3 Action EndBPStudy

End blood pressure study.

#### Syntax

```
public void EndBPStudy()
```

#### Action Errors

|     |                       |
|-----|-----------------------|
| 823 | No SunTech BP monitor |
|-----|-----------------------|

#### 1.1.1.4 Action RecallBPValues

Recall blood pressure values.

#### Syntax

```
public void RecallBPValues()
```

#### Action Errors

|     |                       |
|-----|-----------------------|
| 823 | No SunTech BP monitor |
|-----|-----------------------|

#### 1.1.1.5 Action StartBPSample

Start blood pressure sample.



## Syntax

```
public void StartBPSample()
```

## Action Errors

|     |                       |
|-----|-----------------------|
| 823 | No SunTech BP monitor |
|-----|-----------------------|

### 1.1.1.6 Action StopBPSample

Stop blood pressure sample.

## Syntax

```
public void StopBPSample()
```

## Action Errors

|     |                       |
|-----|-----------------------|
| 823 | No SunTech BP monitor |
|-----|-----------------------|

## 1.2 Service class TreadmillProfileEditor

With the profile editor all stored profiles can be read and the user profiles can be edited and saved. There are three types of profiles in the treadmill:

- Scalable profiles, which can be started in the mode "profile". This profiles are read-only.

Profile numbers:

1 to 6

- Test profiles, which provide specific tests. This profiles are read-only.

Profile numbers:

4 = Bruce

5 = Naughton

6 = Balke

8 = Ellestad A

9 = Ellestad B

21 = User Test Profile 1

22 = User Test Profile 1

23 = User Test Profile 1

24 = User Test Profile 1

70 = MBAFC Interval 600m / 14.5 - 17.5 km/h

71 = MBAFC Interval 500m / 15 - 19 km/h

72 = MBAFC Interval 400m / 16 - 21 km/h

73 = MBAFC Interval 300-600m / 15 - 21 km/h

74 = MBAFC Interval 300-600m / 15 - 21 km/h

75 = MBAFC Interval 300m / 20 km/h

76 = walk protocol 13 steps

80 = VO2 / 10K

81 = VO2 / 11K

82 = VO2 / 12K

84 = VO2 / 14K

85 = SUPER BALKE

90 = Test 90

91 = Test 91

92 = Test 92

93 = Test 93

94 = Test 94

99 = Burn-in

- User profiles, which can be edited and saved.

Profile numbers:

21 to 24 MCU4 (MCU5 28)

Namespace: XDOP.schemas\_coscom\_org.service

Assembly: TreadmillProfileEditor.dll

### Syntax

```
public class TreadmillProfileEditor : IServiceObject
```

### Constructors

|  |  |
|--|--|
| TreadmillProfileEditor(XDOPLib.IService service) | Initializes a new instance of the class. |
|--|--|

### Common Service Properties

| Property name | Value   |
|---------------|---|
| ServiceId     | urn:schemas-coscom-org:serviceId:Profiles               |
| MajorVersion  | 1   |
| MinorVersion  | 1   |
| ServiceType   | urn:schemas-coscom-org:service:TreadmillProfileEditor:1 |

## Service Properties and Events

No Variables

## Service Methods

| Name                                 | Description                                |
|--------------------------------------|--|
| <a href="#">GetProfileHeader</a>     | Gets profile header information.           |
| <a href="#">ReadProfileStep</a>      | Reads profile step information.            |
| <a href="#">WriteUserProfileStep</a> | Writes user test profile step information. |

## Service error codes

| Error code | Error description    |
|------------|----------------------|
| 800        | Not in select mode   |
| 804        | Invalid profile      |
| 813        | Invalid test number  |
| 821        | Invalid step         |
| 822        | Invalid profile type |

## Details

On the following pages you will find a detailed description for every action and variable of this service.

## 1.2.1 Actions

### 1.2.1.1 Action GetProfileHeader

Gets profile header information.

#### Syntax

```
public void GetProfileHeader(  
    System.Byte profileType  
    System.Byte number  
    ref System.Byte steps  
    ref System.Byte intervalType  
    ref System.UInt16 intervalSum  
    ref System.Single maxSpeed  
    ref System.Single maxGrade  
)
```

#### Parameters

##### *profileType*

Profile type:

0 = Factory profile (scalable)

1 = Test profile

2 = User profile

##### *number*

Profile identification number

##### *steps*

Number of steps (max. step index)

##### *intervalType*

Interval type:

0 = mixed

1 = Time

2 = Distance

##### *intervalSum*

Sum of all step intervals (interval type is time or distance)

##### *maxSpeed*

Maximum speed of all steps

##### *maxGrade*

Maximum grade of all steps

#### Remarks

The returned header informations is type-dependent:

Profile type: 0 = Factory profile (scalable)

steps = number of stored steps

intervalType = Time (1) or Distance (2)

maxSpeed = maximum speed of all steps  
maxGrade = maximum grade of all steps

Profile type: 1 = Test profile

steps = number of stored steps  
intervalType = Mixed (0)  
maxSpeed = 0 (no header information available)  
maxGrade = 0 (no header information available)

Profile type: 2 = User profile

steps = 10 (maximum possible number, not the actual number!)  
intervalType = Mixed (0)  
maxSpeed = 0 (no header information available)  
maxGrade = 0 (no header information available)

## Action Errors

|     |                      |
|-----|----------------------|
| 822 | Invalid profile type |
| 804 | Invalid profile      |
| 800 | Not in select mode   |
| 813 | Invalid test number  |

### 1.2.1.2 Action ReadProfileStep

Reads profile step information.

## Syntax

```
public void ReadProfileStep(  
    System.Byte profileType  
    System.Byte number  
    System.Byte step  
    ref System.Byte intervalType  
    ref System.UInt16 interval  
    ref System.Single speed  
    ref System.Byte accelIndex  
    ref System.Single grade  
)
```

### Parameters

#### *profileType*

Profile type:

0 = Factory profile (scalable)  
1 = Test profile  
2 = User profile

#### *number*

Profile identification number

#### *step*

Step index (1 ... max. steps)

#### *intervalType*

Interval type:

1 = Time  
2 = Distance

*interval*  
Time [s] or distance [m]

*speed*  
Speed [m/s]

*accelIndex*  
Acceleration index (1 to 4)

*grade*  
Grade [%]

## Action Errors

|     |                      |
|-----|----------------------|
| 822 | Invalid profile type |
| 804 | Invalid profile      |
| 821 | Invalid step         |
| 800 | Not in select mode   |
| 813 | Invalid test number  |

### 1.2.1.3 Action WriteUserProfileStep

Writes user test profile step information.

## Syntax

```
public void WriteUserProfileStep(  
    System.Byte number  
    System.Byte step  
    System.Byte intervalType  
    System.UInt16 interval  
    System.Single speed  
    System.Byte accelIndex  
    System.Single grade  
)
```

### Parameters

*number*  
Profile identification number

*step*  
Step index (1 ... max. steps)

*intervalType*  
Interval type:

1 = Time  
2 = Distance

*interval*  
Time [s] or distance [m]

*speed*

Speed [m/s]

*accelIndex*

Acceleration index (1 to 4)

*grade*

Grade [%]

## Action Errors

|     |                    |
|-----|--------------------|
| 804 | Invalid profile    |
| 821 | Invalid step       |
| 800 | Not in select mode |

### 1.3 Service class TreadmillGeneralConfig

This service can be used to get general information about the treadmill and to make general configurations.

Namespace: XDOP.schemas\_coscom\_org.service

Assembly: TreadmillGeneralConfig.dll

#### Syntax

```
public class TreadmillGeneralConfig : IServiceObject
```

#### Constructors

|  |  |
|--|--|
| TreadmillGeneralConfig(XDOPLib.IService service) | Initializes a new instance of the class. |
|--|--|

#### Common Service Properties

| Property name | Value   |
|---------------|---|
| ServiceId     | urn:schemas-coscom-org:serviceId:GeneralConfig          |
| MajorVersion  | 1   |
| MinorVersion  | 1   |
| ServiceType   | urn:schemas-coscom-org:service:TreadmillGeneralConfig:1 |

#### Service Properties and Events

| Type        | Name  | Flags | Description   |
|-------------|-------|-------|---|
| System.Byte | Error | E G   | General control error. See action GetErrors for more details. |

#### Service Methods

| Name                               | Description  |
|------------------------------------|--|
| <a href="#">GetDefaultType</a>     | Gets the default type and sub type, which has been set by the factory. This values are important, if the configuration parameters (options) of the treadmill are lost (e.g. the battery of the Real-Time-Clock has been replaced or the signature of the configuration memory is false). |
| <a href="#">GetErrors</a>          | Gets all three types of error codes: user control error, drive error and elevator error.   |
| <a href="#">GetFirmwareVersion</a> | Gets the firmware version.   |
| <a href="#">GetLockMask</a>        | Gets the lock bit mask. Each bit presents a lock for a control mode. If a control mode is locked, the user cannot select and start this mode.  |
| <a href="#">GetOEMCode</a>         | Gets the OEM code of this treadmill. Please contact h/p/cosmos for further details.  |
| <a href="#">GetRTCTime</a>         | Gets the time fields of the battery-buffered real time clock.  |
| <a href="#">GetSerialNumber</a>    | Get the serial number of the treadmill. Contact h/p/cosmos for details.  |



|                                    |  |
|------------------------------------|--|
| <a href="#">GetType</a>            | Gets the type and sub type (treadmill type and variant). Contact h/p/cosmos for a list of all available types and sub types.       |
| <a href="#">SetFailSafeTimeout</a> | This is a very important action, which should be used to set a failsafe timeout before the treadmill is remote controlled.         |
| <a href="#">SetProtocolParams</a>  | Sets optional protocol parameters.   |
| <a href="#">SetRTCTime</a>         | Gets the time fields of the battery-buffered real time clock.  |
| <a href="#">SetSerialNumber</a>    | Get the serial number of the treadmill. Contact h/p/cosmos for details. This action needs two access keys for a successful change. |

## Service error codes

| Error code | Error description |
|------------|-------------------|
| 801        | Wrong access keys |

## Details

On the following pages you will find a detailed description for every action and variable of this service.

### 1.3.1 Actions

#### 1.3.1.1 Action GetDefaultType

Gets the default type and sub type, which has been set by the factory. This values are important, if the configuration parameters (options) of the treadmill are lost (e.g. the battery of the Real-Time-Clock has been replaced or the signature of the configuration memory is false).

#### Syntax

```
public void GetDefaultType(  
    ref System.Byte type  
    ref System.Byte subType  
)
```

Parameters

*type*

Major type (major treadmill type)

*subType*

Sub type (treadmill type variant)

#### 1.3.1.2 Action GetErrors

Gets all three types of error codes: user control error, drive error and elevator error.

#### Syntax

```
public void GetErrors(  
    ref System.Byte controlError  
    ref System.Byte driveError  
    ref System.Byte elevatorError  
)
```

Parameters

*controlError*

General control error

*driveError*

Drive error

*elevatorError*

Elevator error

#### Remarks

Control errors:

0 = No error

1 = Oil interval passed

2 = Service interval passed

20 = Elevator top crash  
21 = Elevator increment error  
30 = Drive increment error  
41 = RTC initialized (backup battery weak)  
50 = FU error

Drive errors:

0 = No error  
1 = Drive increment error (= &gt; control error 30)  
2 = FU error (= &gt; control error 50)

Elevator errors:

0 = No error  
1 = Elevator increment error (= &gt; control error 21)  
2 = Elevator top crash (= &gt; control error 20)

#### 1.3.1.3 Action GetFirmwareVersion

Gets the firmware version.

##### Syntax

```
public void GetFirmwareVersion(  
    ref System.String version  
)
```

Parameters

*version*

A string which is structured like this:  
Major version.subversion 1.subversion 2.Buildnumber.optional appendix

Examples: 4.04.1.0011, 4.04.1.0012.Test

#### 1.3.1.4 Action GetLockMask

Gets the lock bit mask. Each bit presents a lock for a control mode. If a control mode is locked, the user cannot select and start this mode.

##### Syntax

```
public System.Byte[] GetLockMask(  
    ref System.Byte[] mask  
)
```

#### Parameters

*mask*

Lock mask:

0x0F = General

0x10 = Manual

0x20 = Profile

0x40 = Cardio

0x80 = Test

#### 1.3.1.5 Action GetOEMCode

Gets the OEM code of this treadmill. Please contact h/p/cosmos for further details.

#### Syntax

```
public System.Byte GetOEMCode(  
    ref System.Byte code  
)
```

#### Parameters

*code*

OEM code number

#### 1.3.1.6 Action GetRTCTime

Gets the time fields of the battery-buffered real time clock.

#### Syntax

```
public void GetRTCTime(  
    ref System.Byte sec  
    ref System.Byte min  
    ref System.Byte hrs  
    ref System.Byte day  
    ref System.Byte month  
    ref System.UInt16 year  
)
```

#### Parameters

*sec*

RTC seconds: 0 - 59

*min*

RTC minutes: 0 - 59

*hrs*

RTC hours: 0 - 23

*day*

RTC day: 1 - 31

*month*

RTC month: 1 - 12

*year*

RTC year: 2000 - 2099

#### 1.3.1.7 Action GetSerialNumber

Get the serial number of the treadmill. Contact h/p/cosmos for details.

##### Syntax

```
public void GetSerialNumber(  
    ref System.Char code1  
    ref System.Char code2  
    ref System.Char code3  
    ref System.UInt32 num8  
)
```

Parameters

*code1*

Code character 1

*code2*

Code character 2

*code3*

Code character 2

*num8*

Serial number between 0 and 999999999

#### 1.3.1.8 Action GetType

Gets the type and sub type (treadmill type and variant). Contact h/p/cosmos for a list of all available types and sub types.

##### Syntax

```
public void GetType(  
    ref System.Byte type  
    ref System.Byte subType  
)
```

Parameters

*type*

Major type (major treadmill type)

*subType*

Sub type (treadmill type variant)

#### 1.3.1.9 Action SetFailSafeTimeout

This is a very important action, which should be used to set a failsafe timeout before the treadmill is remote controlled.

## Syntax

```
public void SetFailSafeTimeout(  
    System.Byte newTimeout  
    ref System.Byte oldTimeout  
)
```

### Parameters

*newTimeout*

New timeout value in 1/10 seconds (0 to 25.5 seconds)

*oldTimeout*

Old timeout value in 1/10 seconds (0 to 25.5 seconds)

## Remarks

A timeout value of 0 disables the timeout timer and a timeout value between 1 and 255 (0.1 to 25.5 seconds) enables the timer.

The timeout timer increments every 100 milliseconds.

If the set timeout value is reached, the treadmill is stopped (drive and elevator are stopped and the current mode is leaved).

The timeout timer resets and doesn't count up, if the drive is stopped and the elevator doesn't go up or down.

For safety reasons the timeout value should always be set, if the speed or the elevation of the treadmill is changed by an action message.

A reasonable value for the timeout is 1 second (timeout value = 10).

### 1.3.1.10 Action SetProtocolParams

Sets optional protocol parameters.

## Syntax

```
public void SetProtocolParams(  
    System.Boolean processChecksum  
    System.Boolean appendCRLF  
    System.Boolean outputDescNames  
)
```

### Parameters

*processChecksum*

Enables the processing of an 2 hex digit (modulo 256) checksum with element Y0.

*appendCRLF*

Appends the characters CR and LF at the end of a message (after \*Z).

*outputDescNames*

Enables the output of the optional index names in the description lines.

#### 1.3.1.11 Action SetRTCTime

Gets the time fields of the battery-buffered real time clock.

#### Syntax

```
public void SetRTCTime(  
    System.Byte sec  
    System.Byte min  
    System.Byte hrs  
    System.Byte day  
    System.Byte month  
    System.UInt16 year  
)
```

Parameters

*sec*

RTC seconds: 0 - 59

*min*

RTC minutes: 0 - 59

*hrs*

RTC hours: 0 - 23

*day*

RTC day: 1 - 31

*month*

RTC month: 1 - 12

*year*

RTC year: 2000 - 2099

#### 1.3.1.12 Action SetSerialNumber

Get the serial number of the treadmill. Contact h/p/cosmos for details. This action needs two access keys for a successful change.

#### Syntax

```
public void SetSerialNumber(  
    System.Char code1  
    System.Char code2  
    System.Char code3  
    System.UInt32 num8
```

```
System.Byte key1
System.Byte key2
)
```

#### Parameters

*code1*

Code character 1

*code2*

Code character 2

*code3*

Code character 3

*num8*

Serial number between 0 and 99999999

*key1*

Access key 1

*key2*

Access key 2

#### Action Errors

|     |                   |
|-----|-------------------|
| 801 | Wrong access keys |
|-----|-------------------|



## 1.4 Service class TreadmillUserController

This service instantiates a higher level control interface in contrast to the terminal or direct control interface. It models mainly the user control concept of the user terminal. For example one can start a specific test, which is implemented in the treadmill computer.

Namespace: XDOP.schemas\_coscom\_org.service

Assembly: TreadmillUserController.dll

### Syntax

```
public class TreadmillUserController : IServiceObject
```

### Constructors

|   |  |
|---|--|
| TreadmillUserController(XDOPLib.IService service) | Initializes a new instance of the class. |
|---|--|

### Common Service Properties

| Property name | Value  |
|---------------|--|
| ServiceId     | urn:schemas-coscom-org:serviceId: UserControl              |
| MajorVersion  | 1  |
| MinorVersion  | 1  |
| ServiceType   | urn:schemas-coscom-org:service: TreadmillUserController: 1 |

### Service Properties and Events

| Type        | Name              | Flags | Description  |
|-------------|-------------------|-------|--|
| System.Byte | ControlStatus     | E G   | Control status:<br>0 = Stop<br>1 = Run<br>2 = Pause<br>8 = Graded test pause freezed<br>9 = Emergency stop   |
| System.Byte | HrcStatus         | E G   | Actual heart rate control status:<br>0 = idle<br>1 = heart rate too low - load is increased<br>2 = heart rate too high - load is decreased<br>3 = hearte rate is okay<br>4 = failure |
| System.Byte | Mode              | E G   | Actual user control mode:<br>1 = Select<br>2 = Manual<br>3 = Profile<br>4 = Cardio<br>5 = Test<br>6 = User Options<br>7 = Administrator Options                                      |
| System.Byte | NextStepCountdown | E G   | Countdown time value [1/10 s] until next step starts   |
| System.Byte | Step              | E G   | Current step (1 ... 255), if a profile or test is processed  |

### Service Methods

| Name                                   | Description   |
|--|---|
| <a href="#">Cooldown</a>               | Starts the cooldown phase. After the cooldown time (duration) the treadmill stops.  |
| <a href="#">Countdown</a>              | Switches the time display mode from incrementing to decrementing (countdown). The time arguments set the time until the countdown expires. Then the treadmill optionally stops or the time increments again (starts with the the total time). |
| <a href="#">DecGrade</a>               | Decreases the treadmill grade.  |
| <a href="#">DecSpeed</a>               | Decreases the treadmill speed (deceleration).   |
| <a href="#">FreezeStepperTestPause</a> | Freezes the pause of the stepper test (sets the pause time to infinite).  |
| <a href="#">GetHRCStatus</a>           | Gets the actual status of the heart rate control.   |
| <a href="#">GetNextStepStatus</a>      | Gets the status values of the next step. See action ShowNextStep.   |
| <a href="#">GetPersonData</a>          | Gets all person data.   |
| <a href="#">GetProfileStatus</a>       | Gets the actual profile status values.  |
| <a href="#">GetProfListEntry</a>       | Gets a profile entry.   |
| <a href="#">GetProfListSize</a>        | Gets number of stored profiles (profile list size). To get an entry use action GetProfileListEntry.   |
| <a href="#">GetTestListEntry</a>       | Gets a test entry.  |
| <a href="#">GetTestListSize</a>        | Gets number of implemented test (test list size). To get an entry use action GetTestListEntry.  |
| <a href="#">GetTestStatus</a>          | Gets the actual test status, if a test is running.  |
| <a href="#">GetUKKTestResult</a>       | Gets the UKK test result at the end of the test.  |
| <a href="#">HoldGrade</a>              | Holds the actual grade (switches the elevator off).   |
| <a href="#">HoldSpeed</a>              | Holds the actual speed (no acceleration or deceleration).   |
| <a href="#">IncGrade</a>               | Increases the treadmill grade.  |
| <a href="#">IncSpeed</a>               | Increases the treadmill speed (acceleration).   |
| <a href="#">Pause</a>                  | Same as DecSpeed. If the treadmill reaches the speed zero, the speed display shows "PAUS" and the incrementing of the run time is stopped until the treadmill accelerates again.  |
| <a href="#">ProfileStepBack</a>        | Steps one profile step back, if a profile is currently running (profile mode or test mode with a test profile).   |
| <a href="#">ProfileStepForward</a>     | Steps one profile step forward, if a profile is currently running (profile mode or test mode with a test profile).  |
| <a href="#">QuickStop</a>              | Same as action "Stop", but with maximum deceleration.   |
| <a href="#">SetPersonData</a>          | Sets all person data.   |
| <a href="#">SetProfileStepInterval</a> | Sets the interval end time or distance of the current profile step.   |

|  |  |
|--|--|
| <a href="#">ShowNextStep</a>                 | Shows next speed and/or grade values as long as the countdown runs or infinite. At the end of the countdown the actual speed and/or grade are not set. This have to be done with another action. |
| <a href="#">StartHeartRateControl</a>        | Starts heart rate control, which increases or decreases the speed and elevation to get the heart rate between the allowed range (minimum and maximum heart rate).                                |
| <a href="#">StartManual</a>                  | Start manual mode.   |
| <a href="#">StartNextStepperTestInterval</a> | Starts next stepper test intervall (Finishes the pause interval).  |
| <a href="#">StartProfile</a>                 | Start a profile with scale parameters.   |
| <a href="#">StartTestConconi</a>             | Starts the conconi test.   |
| <a href="#">StartTestCooper</a>              | Starts the cooper test.  |
| <a href="#">StartTestGardner</a>             | Starts the gardner test.   |
| <a href="#">StartTestProfile</a>             | Starts a profile-based test.   |
| <a href="#">StartTestRamp</a>                | Starts the ramp test.  |
| <a href="#">StartTestStepper</a>             | Starts stepper test.   |
| <a href="#">StartTestUKK</a>                 | Starts the UKK Test.   |
| <a href="#">Stop</a>                         | Stops the treadmill and the actual mode. It's the same as pressing the key "Stop" at the user terminal.  |
| <a href="#">TestCooperStepBack</a>           | Jumps one cooper test step back. Only valid if a cooper test is running.   |
| <a href="#">TestCooperStepForward</a>        | Jumps one cooper test step forward. Only valid if a cooper test is running.  |

## Service error codes

| Error code | Error description           |
|------------|-----------------------------|
| 802        | No cooldown possible        |
| 803        | No manual mode possible     |
| 804        | Invalid profile             |
| 805        | No OEM profile              |
| 806        | No profile mode possible    |
| 807        | No cardio mode possible     |
| 808        | No test mode possible       |
| 809        | Invalid UKK test parameters |
| 810        | No UKK test active          |
| 811        | No stepper test active      |
| 812        | No cooper test active       |
| 813        | Invalid test number         |
| 814        | Invalid person parameter    |
| 815        | Invalid nextstep mode       |
| 816        | No profile started          |
| 817        | HRC not started             |
| 818        | Test not started            |
| 831        | Speed out of range          |
| 832        | Grade out of range          |

## Details

On the following pages you will find a detailed description for every action and variable of this service.

## 1.4.1 Actions

### 1.4.1.1 Action Cooldown

Starts the cooldown phase. After the cooldown time (duration) the treadmill stops.

#### Syntax

```
public void Cooldown(  
    System.Single speed  
    System.Single grade  
    System.UInt16 duration  
)
```

#### Parameters

*speed*

Cooldown speed [m/s].

*grade*

Cooldown grade [%].

*duration*

Cooldown duration [s].

#### Action Errors

|     |                      |
|-----|----------------------|
| 832 | Grade out of range   |
| 802 | No cooldown possible |
| 831 | Speed out of range   |

### 1.4.1.2 Action Countdown

Switches the time display mode from incrementing to decrementing (countdown). The time arguments set the time until the countdown expires. Then the treadmill optionally stops or the time increments again (starts with the the total time).

#### Syntax

```
public void Countdown(  
    System.Byte hrs  
    System.Byte min  
    System.Byte sec  
    System.Boolean stopFlag  
)
```

#### Parameters

*hrs*

Countdown hours

*min*

Countdown minutes

*sec*

Countdown seconds

*stopFlag*

Flag to stop the treadmill, when the countdown is expired.

#### 1.4.1.3 Action DecGrade

Decreases the treadmill grade.

##### Syntax

```
public void DecGrade()
```

#### 1.4.1.4 Action DecSpeed

Decreases the treadmill speed (deceleration).

##### Syntax

```
public void DecSpeed(  
    System.Byte accelIndex  
)
```

Parameters

*accelIndex*

Index (1 - 7) for a deceleration table with 7 values.

1 = lowest deceleration

...

7 = highest decelartion

#### 1.4.1.5 Action FreezeStepperTestPause

Freezes the pause of the stepper test (sets the pause time to infinite).

##### Syntax

```
public void FreezeStepperTestPause(  
    ref System.Boolean freezed  
)
```

Parameters

*freezed*

#### Action Errors

|     |                       |
|-----|-----------------------|
| 808 | No test mode possible |
|-----|-----------------------|

#### 1.4.1.6 Action GetHRCStatus

Gets the actual status of the heart rate control.

##### Syntax

```

public void GetHRCStatus(
    ref System.Byte heartRate
    ref System.Byte controlState
    ref System.Byte minHeartRate
    ref System.Byte maxHeartRate
)

```

#### Parameters

*heartRate*

Actual heart rate

*controlState*

Actual heart rate control status:

0 = idle

1 = heart rate too low - load is increased

2 = heart rate too high - load is decreased

3 = heart rate is okay

4 = failure

*minHeartRate*

Actual minimum heart rate

*maxHeartRate*

Actual maximum heart rate

## Action Errors

|     |                 |
|-----|-----------------|
| 817 | HRC not started |
|-----|-----------------|

### 1.4.1.7 Action GetNextStepStatus

Gets the status values of the next step. See action ShowNextStep.

## Syntax

```

public void GetNextStepStatus(
    ref System.Byte mode
    ref System.Single speed
    ref System.Single grade
    ref System.Byte countdown
)

```

#### Parameters

*mode*

0 = off

1 = display speed value

2 = display grade value

3 = display both values

*speed*

Speed value [m/s]

*grade*

Grade value [%]

*countdown*

Display time in seconds.

0 = show next values infinite

1 - 25 = with countdown peeps and automatic clearing of the shown values

#### 1.4.1.8 Action GetPersonData

Gets all person data.

##### Syntax

```
public void GetPersonData(  
    ref System.Char sex  
    ref System.Byte age  
    ref System.Byte weight  
    ref System.Byte size  
)
```

Parameters

*sex*

Person sex:

'F' = female

'M' = make

*age*

Person age [years]

*weight*

Person weight [kg]

*size*

Person size [cm]

#### 1.4.1.9 Action GetProfileStatus

Gets the actual profile status values.

##### Syntax

```
public void GetProfileStatus(  
    ref System.Byte number  
    ref System.Byte step  
    ref System.Single speed  
    ref System.Single grade  
    ref System.Byte intervalType  
    ref System.UInt16 intervalStart  
    ref System.UInt16 intervalEnd  
)
```

Parameters

*number*

Identification number of the profile



*step*

Current step number (1 ... 255)

*speed*

Actual speed [m/s]

*grade*

Actual grade [%]

*intervalType*

Interval type:

1 = time

2 = distance

*intervalStart*

Interval start time [s] or distance [m]

*intervalEnd*

Interval end time [s] or distance [m]

## Action Errors

|     |                    |
|-----|--------------------|
| 816 | No profile started |
|-----|--------------------|

### 1.4.1.10 Action GetProfListEntry

Gets a profile entry.

## Syntax

```
public void GetProfListEntry(  
    System.Byte index  
    ref System.Byte number  
    ref System.Byte intervalType  
    ref System.Boolean allowed  
)
```

Parameters

*index*

Profile list index (0 to GetProfileListSize() - 1)

*number*

Displayed (identification) number

*intervalType*

Profile interval type:

0 = time and distance steps

1 = time steps

2 = distance steps

*allowed*

True, if this profile can be started

False, if this profile is not allowed by the actual treadmill configuration

#### 1.4.1.11 Action GetProfListSize

Gets number of stored profiles (profile list size). To get an entry use action GetProfileListEntry.

##### Syntax

```
public System.Byte GetProfListSize(  
    ref System.Byte size  
)
```

##### Parameters

*size*

Number of stored profiles

#### 1.4.1.12 Action GetTestListEntry

Gets a test entry.

##### Syntax

```
public void GetTestListEntry(  
    System.Byte index  
    ref System.String name  
    ref System.Byte number  
    ref System.Byte testType  
    ref System.Boolean allowed  
)
```

##### Parameters

*index*

Test list index (0 to GetTestListSize() - 1)

*name*

Test name

*number*

Displayed (identification) number

*testType*

Test type:

1 = complex

2 = factory profile

3 = user profile

*allowed*

True, if this test can be started

False, if this test is not allowed by the actual treadmill configuration

#### 1.4.1.13 Action GetTestListSize

Gets number of implemented test (test list size). To get an entry use action GetTestListEntry.

## Syntax

```
public System.Byte GetTestListSize(  
    ref System.Byte size  
)
```

### Parameters

*size*

Number of implemented tests

#### 1.4.1.14 Action GetTestStatus

Gets the actual test status, if a test is running.

## Syntax

```
public void GetTestStatus(  
    ref System.String name  
    ref System.Byte number  
    ref System.Byte testType  
    ref System.Byte step  
)
```

### Parameters

*name*

Test name

*number*

Test identification number

*testType*

Test type:

1 = Complex

2 = Profile

3 = User Profile

*step*

Current step number (1 ... 255)

## Action Errors

|     |                  |
|-----|------------------|
| 818 | Test not started |
|-----|------------------|

#### 1.4.1.15 Action GetUKKTestResult

Gets the UKK test result at the end of the test.

## Syntax

```

public void GetUKKTestResult(
    ref System.Byte heartRate500
    ref System.Byte heartRate1000
    ref System.Byte heartRate1500
    ref System.Byte heartRate2000
    ref System.Byte heartRateAverage
    ref System.UInt16 time2000
    ref System.Byte index
)

```

#### Parameters

*heartRate500*

Heart Rate at 500 meters

*heartRate1000*

Heart Rate at 1000 meters

*heartRate1500*

Heart Rate at 1500 meters

*heartRate2000*

Heart Rate at 2000 meters

*heartRateAverage*

Heart rate average

*time2000*

Time [s] for test length of 2000 meters

*index*

UKK test index

Interpretation of the index:

< 70: Considerably below average fitness

70 - 89: Slightly below average fitness

90 - 110: Average fitness

111 - 130: Slightly above average fitness

> 130: Considerably above average fitness

## Action Errors

|     |                    |
|-----|--------------------|
| 810 | No UKK test active |
|-----|--------------------|

### 1.4.1.16 Action HoldGrade

Holds the actual grade (switches the elevator off).

## Syntax

```

public void HoldGrade()

```

### 1.4.1.17 Action HoldSpeed

Holds the actual speed (no acceleration or deceleration).

## Syntax

```
public void HoldSpeed()
```

### 1.4.1.18 Action IncGrade

Increases the treadmill grade.

## Syntax

```
public void IncGrade()
```

### 1.4.1.19 Action IncSpeed

Increases the treadmill speed (acceleration).

## Syntax

```
public void IncSpeed(  
    System.Byte accelIndex  
)
```

Parameters

*accelIndex*

Index (1 - 7) for an acceleration table with 7 values.

1 = lowest acceleration

...

7 = highest acceleration

### 1.4.1.20 Action Pause

Same as DecSpeed. If the treadmill reaches the speed zero, the speed display shows "PAUS" and the incrementing of the run time is stopped until the treadmill accelerates again.

## Syntax

```
public void Pause()
```

### 1.4.1.21 Action ProfileStepBack

Steps one profile step back, if a profile is currently running (profile mode or test mode with a test profile).

## Syntax

```
public void ProfileStepBack()
```

## Action Errors

|     |                          |
|-----|--------------------------|
| 806 | No profile mode possible |
|-----|--------------------------|

### 1.4.1.22 Action ProfileStepForward

Steps one profile step forward, if a profile is currently running (profile mode or test mode with a test profile).

## Syntax

```
public void ProfileStepForward()
```

## Action Errors

|     |                          |
|-----|--------------------------|
| 806 | No profile mode possible |
|-----|--------------------------|

### 1.4.1.23 Action QuickStop

Same as action "Stop", but with maximum deceleration.

## Syntax

```
public void QuickStop()
```

### 1.4.1.24 Action SetPersonData

Sets all person data.

## Syntax

```
public void SetPersonData(  
    System.Char sex  
    System.Byte age  
    System.Byte weight  
    System.Byte size  
)
```

### Parameters

*sex*

Person sex:

'F' = female

'M' = male

*age*

Person age [years]

*weight*

Person weight [kg]

*size*

Person size [cm]

## Action Errors

|     |                          |
|-----|--------------------------|
| 814 | Invalid person parameter |
|-----|--------------------------|

### 1.4.1.25 Action SetProfileStepInterval

Sets the interval end time or distance of the current profile step.

## Syntax

```
public void SetProfileStepInterval(  
    System.UInt16 intervalEnd  
)
```

Parameters

*intervalEnd*

Interval end time [s] or distance [m]

## Action Errors

|     |                    |
|-----|--------------------|
| 816 | No profile started |
|-----|--------------------|

### 1.4.1.26 Action ShowNextStep

Shows next speed and/or grade values as long as the countdown runs or infinite. At the end of the countdown the actual speed and/or grade are not set. This have to be done with another action.

## Syntax

```
public void ShowNextStep(  
    System.Byte mode  
    System.Single speed  
    System.Single grade  
    System.Byte time  
)
```

Parameters

*mode*

The mode determines which values are shown:

0 = off (no argument values are shown and a running countdown is canceled)

1 = display speed value

2 = display grade value

3 = display both values

*speed*

Speed value [m/s]

*grade*

Grade value [%]

*time*

Display time in seconds.

0 = show next values infinite

1 - 25 = with countdown peeps and automatic clearing of the shown values

## Action Errors

|     |                       |
|-----|-----------------------|
| 815 | Invalid nextstep mode |
|-----|-----------------------|

### 1.4.1.27 Action StartHeartRateControl

Starts heart rate control, which increases or decreases the speed and elevation to get the heart rate between the allowed range (minimum and maximum heart rate).

## Syntax

```
public void StartHeartRateControl(  
    System.Single maxSpeed  
    System.Byte maxHeartrate  
    System.Byte minHeartrate  
)
```

### Parameters

#### *maxSpeed*

Maximum speed to get a higher heart rate. If this speed value isn't sufficient for the heart rate control, then the grade is increased.

#### *maxHeartrate*

Maximum heart rate

#### *minHeartrate*

Minimum heart rate

## Action Errors

|     |                         |
|-----|-------------------------|
| 831 | Speed out of range      |
| 807 | No cardio mode possible |

### 1.4.1.28 Action StartManual

Start manual mode.

## Syntax

```
public void StartManual()
```

## Action Errors

|     |                         |
|-----|-------------------------|
| 803 | No manual mode possible |
|-----|-------------------------|

### 1.4.1.29 Action StartNextStepperTestInterval

Starts next stepper test intervall (Finishes the pause interval).

## Syntax



```
public void StartNextStepperTestInterval()
```

## Action Errors

|     |                        |
|-----|------------------------|
| 811 | No stepper test active |
|-----|------------------------|

### 1.4.1.30 Action StartProfile

Start a profile with scale parameters.

## Syntax

```
public void StartProfile(  
    System.Byte number  
    System.Byte scaleSpeed  
    System.Byte scaleGrade  
    System.Byte scaleInterval  
)
```

### Parameters

#### *number*

Profile identification (displayed) number

#### *scaleSpeed*

Speed scale parameter. Allowd values:

60 = 60 %

80 = 80 %

100 = 100 %

120 = 120 %

140 = 140 %

160 = 160 %

#### *scaleGrade*

Grade scale parameter. Allowd values:

60 = 60 %

80 = 80 %

100 = 100 %

120 = 120 %

140 = 140 %

160 = 160 %

#### *scaleInterval*

Interval scale parameter. Allowd values:

60 = 60 %

80 = 80 %

100 = 100 %

120 = 120 %

140 = 140 %

160 = 160 %

## Action Errors

|     |                          |
|-----|--------------------------|
| 804 | Invalid profile          |
| 806 | No profile mode possible |
| 805 | No OEM profile           |

#### 1.4.1.31 Action StartTestConconi

Starts the conconi test.

#### Syntax

```
public void StartTestConconi(
    System.Single startSpeed
    System.UInt16 interval
    System.Single speedIncrement
)
```

Parameters

*startSpeed*

Start speed [m/s]

*interval*

Interval distance [m] of test steps

*speedIncrement*

Speed increment [m/s] between the test steps

#### Action Errors

|     |                       |
|-----|-----------------------|
| 808 | No test mode possible |
| 831 | Speed out of range    |

#### 1.4.1.32 Action StartTestCooper

Starts the cooper test.

#### Syntax

```
public void StartTestCooper()
```

#### Action Errors

|     |                       |
|-----|-----------------------|
| 808 | No test mode possible |
|-----|-----------------------|

#### 1.4.1.33 Action StartTestGardner

Starts the gardner test.

#### Syntax

```
public void StartTestGardner()
```

#### Action Errors

|     |                       |
|-----|-----------------------|
| 808 | No test mode possible |
|-----|-----------------------|

#### 1.4.1.34 Action StartTestProfile

Starts a profile-based test.

##### Syntax

```
public void StartTestProfile(  
    System.Byte number  
)
```

Parameters

*number*

Test identification number

##### Action Errors

|     |                       |
|-----|-----------------------|
| 808 | No test mode possible |
| 813 | Invalid test number   |

#### 1.4.1.35 Action StartTestRamp

Starts the ramp test.

##### Syntax

```
public void StartTestRamp(  
    System.Single endSpeed  
    System.Byte interval  
)
```

Parameters

*endSpeed*

End speed [m/s]

*interval*

Interval time [s]

##### Action Errors

|     |                       |
|-----|-----------------------|
| 808 | No test mode possible |
| 831 | Speed out of range    |

#### 1.4.1.36 Action StartTestStepper

Starts stepper test.

##### Syntax

```
public void StartTestStepper(  
    System.Single startSpeed  
    System.UInt16 intervalTime  
    System.Single speedIncrement  
    System.Byte accelIndex  
    System.UInt16 pauseTime  
)
```

#### Parameters

*startSpeed*

Start speed [m/s]

*intervalTime*

Interval time [s]

*speedIncrement*

Speed increment [m/s]

*accelIndex*

Acceleration/Deceleration index (1 to 7)

*pauseTime*

Pause time [s]

#### Action Errors

|     |                       |
|-----|-----------------------|
| 808 | No test mode possible |
| 831 | Speed out of range    |

#### 1.4.1.37 Action StartTestUKK

Starts the UKK Test.

#### Syntax

```
public void StartTestUKK(  
    System.Char sex  
    System.Byte age  
    System.Byte weight  
    System.Byte size  
)
```

#### Parameters

*sex*

Sex of the running person:

'F' = female

'M' = male

*age*

Person age [years]

*weight*

Person weight [kg]

*size*

Person size [cm]

#### Action Errors

|     |                             |
|-----|-----------------------------|
| 808 | No test mode possible       |
| 809 | Invalid UKK test parameters |

#### 1.4.1.38 Action Stop

Stops the treadmill and the actual mode. It's the same as pressing the key "Stop" at the user terminal.

##### Syntax

```
public void Stop()
```

#### 1.4.1.39 Action TestCooperStepBack

Jumps one cooper test step back. Only valid if a cooper test is running.

##### Syntax

```
public void TestCooperStepBack()
```

##### Action Errors

|     |                       |
|-----|-----------------------|
| 812 | No cooper test active |
|-----|-----------------------|

#### 1.4.1.40 Action TestCooperStepForward

Jumps one cooper test step forward. Only valid if a cooper test is running.

##### Syntax

```
public void TestCooperStepForward()
```

##### Action Errors

|     |                       |
|-----|-----------------------|
| 812 | No cooper test active |
|-----|-----------------------|

## 1.5 Service class TreadmillOptions

This is a management service for all optional configuration settings and values.

Namespace: XDOP.schemas\_coscom\_org.service

Assembly: TreadmillOptions.dll

### Syntax

```
public class TreadmillOptions : IServiceObject
```

### Constructors

|   |  |
|---|--|
| <code>TreadmillOptions(XDOPLib.IService service)</code> | Initializes a new instance of the class. |
|---|--|

### Common Service Properties

| Property name | Value  |
|---------------|--|
| ServiceId     | urn:schemas-coscom-org:serviceId:Options           |
| MajorVersion  | 1  |
| MinorVersion  | 1  |
| ServiceType   | urn:schemas-coscom-org:service:TreadmillOptions: 1 |

### Service Properties and Events

No Variables

### Service Methods

| Name                         | Description                           |
|------------------------------|---------------------------------------|
| <a href="#">GetListEntry</a> | Gets an entry from a list.            |
| <a href="#">GetListSize</a>  | Gets the number of entries in a list. |
| <a href="#">GetOption</a>    | Gets an option value and its range.   |
| <a href="#">GetVersion</a>   | Gets the version of the options.      |
| <a href="#">SetOption</a>    | Sets an option value.                 |

### Service error codes

| Error code | Error description                                |
|------------|--|
| 800        | Not in select mode                               |
| 801        | Wrong access keys                                |
| 824        | Unknown option number                            |
| 825        | Unknown option error                             |
| 826        | Option value out of range                        |
| 827        | Invalid option list index                        |
| 828        | Invalid option level                             |
| 829        | No option list available                         |
| 833        | Current communication protocol cannot be changed |

### Details

On the following pages you will find a detailed description for every action and variable of this service.

## 1.5.1 Actions

### 1.5.1.1 Action GetListEntry

Gets an entry from a list.

#### Syntax

```
public void GetListEntry(  
    System.Char level  
    System.Byte number  
    System.Byte index  
    ref System.Int32 _value  
)
```

#### Parameters

*level*

Option level:

'A' = Administrator

'U' = User

*number*

Option identification number or 0 for a list with all options

*index*

List index (0 to GetListSize() - 1)

*\_value*

List entry value

#### Action Errors

|     |                           |
|-----|---------------------------|
| 828 | Invalid option level      |
| 824 | Unknown option number     |
| 829 | No option list available  |
| 800 | Not in select mode        |
| 825 | Unknown option error      |
| 827 | Invalid option list index |

### 1.5.1.2 Action GetListSize

Gets the number of entries in a list.

#### Syntax

```
public void GetListSize(  
    System.Char level  
    System.Byte number  
    ref System.Byte size  
)
```

#### Parameters

*level*

Option level:

'A' = Administrator

'U' = User

*number*

Option identification number or 0 for options list

*size*

Number of entries in the list

## Action Errors

|     |                          |
|-----|--------------------------|
| 828 | Invalid option level     |
| 824 | Unknown option number    |
| 829 | No option list available |
| 800 | Not in select mode       |
| 825 | Unknown option error     |

### 1.5.1.3 Action GetOption

Gets an option value and its range.

## Syntax

```
public void GetOption(  
    System.Char level  
    System.Byte number  
    ref System.Int32 _value  
    ref System.Int32 min  
    ref System.Int32 max  
)
```

Parameters

*level*

Option level:

'A' = Administrator

'U' = User

*number*

Option identification number

*\_value*

Actual option value

*min*

Minimum option value

*max*

Maximum option value

## Action Errors

|     |                       |
|-----|-----------------------|
| 828 | Invalid option level  |
| 824 | Unknown option number |
| 800 | Not in select mode    |
| 825 | Unknown option error  |



#### 1.5.1.4 Action GetVersion

Gets the version of the options.

##### Syntax

```
public System.String GetVersion(  
    ref System.String version  
)
```

Parameters

*version*

Version string e.g. "101"

#### 1.5.1.5 Action SetOption

Sets an option value.

##### Syntax

```
public void SetOption(  
    System.Char level  
    System.Byte number  
    System.Int32 _value  
    System.Byte key1  
    System.Byte key2  
)
```

Parameters

*level*

Option level:

'A' = Administrator

'U' = User

*number*

Option identification number

*\_value*

New option value

*key1*

Access key 1

*key2*

Access key 2

##### Action Errors

|     |  |
|-----|--|
| 828 | Invalid option level                             |
| 824 | Unknown option number                            |
| 800 | Not in select mode                               |
| 826 | Option value out of range                        |
| 825 | Unknown option error                             |
| 833 | Current communication protocol cannot be changed |

|     |                   |
|-----|-------------------|
| 801 | Wrong access keys |
|-----|-------------------|

## 1.6 Service class TreadmillDirectControl

The direct control service has actions to control the speed of the drive and the grade of the elevator directly.

Namespace: XDOP.schemas\_coscom\_org.service

Assembly: TreadmillDirectControl.dll

### Syntax

```
public class TreadmillDirectControl : IServiceObject
```

### Constructors

|  |  |
|--|--|
| TreadmillDirectControl(XDOPLib.IService service) | Initializes a new instance of the class. |
|--|--|

### Common Service Properties

| Property name | Value   |
|---------------|---|
| ServiceId     | urn:schemas-coscom-org:serviceId:DirectControl          |
| MajorVersion  | 1   |
| MinorVersion  | 1   |
| ServiceType   | urn:schemas-coscom-org:service:TreadmillDirectControl:1 |

### Service Properties and Events

| Type          | Name           | Flags | Description   |
|---------------|----------------|-------|---|
| System.Single | Acceleration   | E G   | Actual acceleration (>0) deceleration (<0) value [m/s <sup>2</sup> ]  |
| System.Single | ActualGrade    | E G   | Actual grade value [%]  |
| System.Single | ActualSpeed    | E G   | Actual speed value [m/s]  |
| System.Byte   | DriveStatus    | E G   | Actual drive status:<br>0 = Stopped (RFR off)<br>1 = Constant speed<br>2 = Acceleration<br>3 = Deceleration |
| System.Byte   | ElevatorStatus | E G   | Actual elevator status:<br>0 = Off<br>1 = Up<br>2 = Down  |
| System.Single | FinalGrade     | E G   | Final grade value [%]   |
| System.Single | FinalSpeed     | E G   | Final speed value [m/s]   |

### Service Methods

| Name                                    | Description  |
|---|--|
| <a href="#">GetDriveAccelDecelRange</a> | Gets the configured acceleration/deceleration range.           |
| <a href="#">GetDriveSpeedRange</a>      | Gets the configured speed range.                               |
| <a href="#">GetDriveStatus</a>          | Gets drive status values.                                      |
| <a href="#">GetElevatorGradeRange</a>   | Gets the configured grade range.                               |
| <a href="#">GetElevatorStatus</a>       | Gets elevator status values.                                   |
| <a href="#">SetDriveSpeed</a>           | Sets a new drive speed and an acceleration/deceleration value. |
| <a href="#">SetElevatorGrade</a>        | Sets a new grade value.  |

## Service error codes

| Error code | Error description  |
|------------|--------------------|
| 831        | Speed out of range |
| 832        | Grade out of range |

## Details

On the following pages you will find a detailed description for every action and variable of this service.

## 1.6.1 Actions

### 1.6.1.1 Action GetDriveAccelDecelRange

Gets the configured acceleration/deceleration range.

#### Syntax

```
public void GetDriveAccelDecelRange(  
    ref System.Byte minIndex  
    ref System.Byte maxIndex  
    ref System.Single minValue  
    ref System.Single maxValue  
)
```

#### Parameters

*minIndex*

Minimum acceleration/deceleration index value

*maxIndex*

Maximum acceleration/deceleration index value

*minValue*

Minimum acceleration/deceleration value [m/s<sup>2</sup>]

*maxValue*

Maximum acceleration/deceleration value [m/s<sup>2</sup>]

### 1.6.1.2 Action GetDriveSpeedRange

Gets the configured speed range.

#### Syntax

```
public void GetDriveSpeedRange(  
    ref System.Single minSpeed  
    ref System.Single maxSpeed  
)
```

#### Parameters

*minSpeed*

Minimum speed [m/s]

*maxSpeed*

Maximum speed [m/s]

### 1.6.1.3 Action GetDriveStatus

Gets drive status values.

#### Syntax

```
public void GetDriveStatus(  
    ref System.Byte driveStatus  
    ref System.Single actualSpeed
```

```

    ref System.Single finalSpeed
    ref System.Single acceleration
)

```

#### Parameters

##### *driveStatus*

Drive status:

- 0 = Stopped (RFR off)
- 1 = Constant speed
- 2 = Acceleration
- 3 = Deceleration

##### *actualSpeed*

Actual Speed [m/s]

##### *finalSpeed*

Final speed [m/s]

##### *acceleration*

Acceleration [m/s<sup>2</sup>]

### 1.6.1.4 Action GetElevatorGradeRange

Gets the configured grade range.

#### Syntax

```

public void GetElevatorGradeRange(
    ref System.Single minGrade
    ref System.Single maxGrade
)

```

#### Parameters

##### *minGrade*

Minimum grade [%]

##### *maxGrade*

Maximum grade [%]

### 1.6.1.5 Action GetElevatorStatus

Gets elevator status values.

#### Syntax

```

public void GetElevatorStatus(
    ref System.Byte elevatorStatus
    ref System.Single actualGrade
    ref System.Single finalGrade
)

```

#### Parameters

##### *elevatorStatus*

Elevator status:

- 0 = Off

1 = Up  
2 = Down

*actualGrade*  
Actual Grade [%]

*finalGrade*  
Final Grade [%]

#### 1.6.1.6 Action SetDriveSpeed

Sets a new drive speed and an acceleration/deceleration value.

##### Syntax

```
public void SetDriveSpeed(  
    System.Single speed  
    System.Byte accelIndex  
    System.Single acceleration  
)
```

##### Parameters

*speed*  
Final Speed [m/s]

*accelIndex*  
Acceleration index:  
0 = Use "acceleration" argument  
1 = (Max. FU speed) / 131 s  
2 = (Max. FU speed) / 65.5 s  
3 = (Max. FU speed) / 32.8 s  
4 = (Max. FU speed) / 16.0 s  
5 = (Max. FU speed) / 8.0 s  
6 = (Max. FU speed) / 5.0 s  
7 = (Max. FU speed) / 3.0 s

*acceleration*  
Acceleration [m/s<sup>2</sup>]

##### Action Errors

|     |                    |
|-----|--------------------|
| 831 | Speed out of range |
|-----|--------------------|

#### 1.6.1.7 Action SetElevatorGrade

Sets a new grade value.

##### Syntax

```
public void SetElevatorGrade(  
    System.Single grade  
)
```

Parameters

*grade*

Final Grade [%]

## Action Errors

|     |                    |
|-----|--------------------|
| 832 | Grade out of range |
|-----|--------------------|



## 1.7 Service class TerminalLCD6x4Panel

This is a service to build a remote terminal with the same display and keyboard layout as the built-in treadmill terminal.

Namespace: XDOP.schemas\_coscom\_org.service

Assembly: TerminalLCD6x4Panel.dll

### Syntax

```
public class TerminalLCD6x4Panel : IServiceObject
```

### Constructors

|   |  |
|---|--|
| TerminalLCD6x4Panel(XDOPLib.IService service) | Initializes a new instance of the class. |
|---|--|

### Common Service Properties

| Property name | Value  |
|---------------|--|
| ServiceId     | urn:schemas-coscom-org:serviceId:Terminal            |
| MajorVersion  | 1  |
| MinorVersion  | 1  |
| ServiceType   | urn:schemas-coscom-org:service:TerminalLCD6x4Panel:1 |

### Service Properties and Events

| Type           | Name       | Flags | Description   |
|----------------|------------|-------|---|
| System.Boolean | BlinkState | E G   | Blink state of all LCD displays and LEDs. Every change sends an event.<br>0 = Off, 1 = On   |
| System.Byte    | BuzzerTime | E G   | Current remaining buzzer on time in 1/100 seconds. The change events are sent as often as possible.   |
| System.Byte[]  | KeyState   | E G   | Bit field with MCU 4 key states (key on = bit set, key off = bit cleared).<br>Key bit masks:<br>0x01 = STOP<br>0x02 = START<br>0x04 = UP<br>0x08 = DOWN<br>0x10 = PLUS<br>0x20 = MINUS<br><br>0x80 = EMERGENCY STOP<br><br>Each state change sends an event, if the change rate isn't too high. |

|               |           |     |  |
|---------------|-----------|-----|--|
| System.Byte[] | LEDGroups | E G | <p>9 hex digits with the LED index of a LED group.<br/> Indexes of LED groups:<br/> first digit:<br/> 5 = manual<br/> 4 = profile<br/> 3 = cardio<br/> 2 = test</p> <p>2. digit<br/> 5 = max</p> <p>3. digit<br/> 1 = mph<br/> 2 = m/s<br/> 3 = km/h<br/> 4 = m/min</p> <p>4. digit<br/> 3 = miles<br/> 4 = km<br/> 5 = m</p> <p>5. digit<br/> 4 = &lt;°<br/> 5 = %</p> <p>6. digit<br/> 1 = No.<br/> 2 = Step</p> <p>7. digit<br/> 1 = power<br/> 2 = energy<br/> 5 = met</p> <p>8. digit<br/> 4 = v min.<br/> 5 = ^ max.</p> <p>9. digit<br/> 1 = weight<br/> 2 = age<br/> 3 = sex</p> <p>index + A = blinking</p> <p>0 = all off</p> <p>Each change sends an event.</p> |
|---------------|-----------|-----|--|

|               |           |     |   |
|---------------|-----------|-----|---|
| System.String | LLDisplay | E G | <p>Lower left LCD display string and attributes:</p> <ul style="list-style-type: none"> <li>- the first 4 characters build the display string</li> <li>- the fifth character (hexdecimal) describes the bit field for blinking digits</li> <li>- the sixth character (hexdecimal) describes the bit field for points between the digits</li> <li>- the seventh character (hexdecimal) describes the bit field for blinking points</li> </ul> <p>Bit pattern for blinking display digits (character 5 in the message string)</p> <p>0x0 No blinking digits<br/> 0x1 Left digit is blinking<br/> 0x2 Middle left digit is blinking<br/> 0x4 Middle right digit is blinking<br/> 0x8 Right digit is blinking</p> <p>Bit pattern for display points (character 6 in the message string)</p> <p>0x0 No points<br/> 0x1 Left point<br/> 0x2 Middle point<br/> 0x4 Right point<br/> 0x8 Colon</p> <p>Bit pattern for blinking display points (character 7 in the message string)</p> <p>0x0 No blinking points<br/> 0x1 Blinking left point<br/> 0x2 Blinking middle point<br/> 0x4 Blinking right point<br/> 0x8 Blinking colon</p> <p>Example:<br/> " 00040" shows the string " 00" with a right point and looks like " 0.0"</p> <p>Each state change sends an event, if the change rate isn't too high.</p> |
| System.String | LMDisplay | E G | Lower middle LCD display string and attributes. Same format as LLLDisplay.  |
| System.String | LRDisplay | E G | Lower right LCD display string and attributes. Same format as LLLDisplay.   |
| System.String | ULDisplay | E G | Upper left LCD display string and attributes. Same format as LLLDisplay.  |
| System.String | UMDisplay | E G | Upper middle LCD display string and attributes. Same format as LLLDisplay.  |
| System.String | URDisplay | E G | Upper right LCD display string and attributes. Same format as LLLDisplay.   |

## Service Methods

| Name                 | Description  |
|----------------------|--|
| <a href="#">Beep</a> | Switches the buzzer of the treadmill terminal on. The beep lasts as long as the given duration argument. |

|                                 |  |
|---------------------------------|--|
| <a href="#">GetDistanceUnit</a> | Gets the configured distance unit.   |
| <a href="#">GetGradeUnit</a>    | Gets the configured grade unit.  |
| <a href="#">GetKeyLock</a>      | Gets the keyboard lock state. If locked, only the key "-" and "Stop" can be used. Both stop the treadmill.                               |
| <a href="#">GetSpeedUnit</a>    | Gets the configured speed unit.  |
| <a href="#">KeyDown</a>         | Changes the key state to "key down" (key pressed).   |
| <a href="#">KeyUp</a>           | Changes the key state to "key up" (key released).  |
| <a href="#">SetDistanceUnit</a> | Sets the distance unit temporarily. If the treadmill is switched off or if the configuration mode is entered, then this setting is lost. |
| <a href="#">SetGradeUnit</a>    | Sets the grade unit temporarily. If the treadmill is switched off or if the configuration mode is entered, then this setting is lost.    |
| <a href="#">SetKeyLock</a>      | Sets the keyboard lock state. If locked, only the key "-" and "Stop" can be used. Both stop the treadmill.                               |
| <a href="#">SetModeLED</a>      | Sets one of the mode LEDs or switches all off. The current control mode isn't be influenced.   |
| <a href="#">SetSpeedUnit</a>    | Sets the speed unit temporarily. If the treadmill is switched off or if the configuration mode is entered, then this setting is lost.    |

## Service error codes

| Error code | Error description |
|------------|-------------------|
| 819        | Invalid LED mode  |
| 820        | Invalid unit      |

## Details

On the following pages you will find a detailed description for every action and variable of this service.

### 1.7.1 Actions

#### 1.7.1.1 Action Beep

Switches the buzzer of the treadmill terminal on. The beep lasts as long as the given duration argument.

##### Syntax

```
public void Beep(  
    System.Byte duration  
)
```

Parameters

*duration*

Beep time in 1/100 seconds (0 to 2.55 s).

#### 1.7.1.2 Action GetDistanceUnit

Gets the configured distance unit.

##### Syntax

```
public System.Byte GetDistanceUnit(  
    ref System.Byte unit  
)
```

Parameters

*unit*

0 = km, 1 = miles, 2 = m

#### 1.7.1.3 Action GetGradeUnit

Gets the configured grade unit.

##### Syntax

```
public System.Byte GetGradeUnit(  
    ref System.Byte unit  
)
```

Parameters

*unit*

0 = %, 1 = ° (degree)

#### 1.7.1.4 Action GetKeyLock

Gets the keyboard lock state. If locked, only the key "-" and "Stop" can be used. Both stop the treadmill.

##### Syntax

```
public void GetKeyLock(  
    ref System.Boolean _lock  
)
```

Parameters

*\_lock*  
0 = unlocked, 1 = locked

#### 1.7.1.5 Action GetSpeedUnit

Gets the configured speed unit.

##### Syntax

```
public System.Byte GetSpeedUnit(  
    ref System.Byte unit  
)
```

Parameters

*unit*  
0 = km/h, 1 = m/s, 2 = mph, 3 = m/min

#### 1.7.1.6 Action KeyDown

Changes the key state to "key down" (key pressed).

##### Syntax

```
public void KeyDown(  
    System.Char key  
)
```

Parameters

*key*  
+ = Minus, - = Plus, U = Up, D = Down, E = Start, S = Stop, F = Emergency Stop

#### 1.7.1.7 Action KeyUp

Changes the key state to "key up" (key released).

##### Syntax

```
public void KeyUp(  
    System.Char key  
)
```

Parameters

*key*  
+ = Minus, - = Plus, U = Up, D = Down, E = Start, S = Stop, F = Emergency Stop

#### 1.7.1.8 Action SetDistanceUnit

Sets the distance unit temporarily. If the treadmill is switched off or if the configuration mode is entered, then this setting is lost.

##### Syntax

```
public void SetDistanceUnit(  
    System.Byte unit  
)
```

Parameters

*unit*

0 = km, 1 = miles, 2 = m

##### Action Errors

|     |              |
|-----|--------------|
| 820 | Invalid unit |
|-----|--------------|

#### 1.7.1.9 Action SetGradeUnit

Sets the grade unit temporarily. If the treadmill is switched off or if the configuration mode is entered, then this setting is lost.

##### Syntax

```
public void SetGradeUnit(  
    System.Byte unit  
)
```

Parameters

*unit*

0 = %, 1 = ° (degree)

##### Action Errors

|     |              |
|-----|--------------|
| 820 | Invalid unit |
|-----|--------------|

#### 1.7.1.10 Action SetKeyLock

Sets the keyboard lock state. If locked, only the key "-" and "Stop" can be used. Both stop the treadmill.

##### Syntax

```
public void SetKeyLock(  
    System.Boolean _lock  
)
```

Parameters

*\_lock*

0 = unlocked, 1 = locked

#### 1.7.1.11 Action SetModeLED

Sets one of the mode LEDs or switches all off. The current control mode isn't be influenced.

##### Syntax

```
public void SetModeLED(  
    System.Byte mode  
    System.Boolean blink  
)
```

Parameters

*mode*

0 = All Off, 1 = Manual, 2 = Profile, 3 = Cardio, 4 = Test

*blink*

0 = permanent on, 1 = blinking

##### Action Errors

|     |                  |
|-----|------------------|
| 819 | Invalid LED mode |
|-----|------------------|

#### 1.7.1.12 Action SetSpeedUnit

Sets the speed unit temporarily. If the treadmill is switched off or if the configuration mode is entered, then this setting is lost.

##### Syntax

```
public void SetSpeedUnit(  
    System.Byte unit  
)
```

Parameters

*unit*

0 = km/h, 1 = m/s, 2 = mph, 3 = m/min

##### Action Errors

|     |              |
|-----|--------------|
| 820 | Invalid unit |
|-----|--------------|



## 1.8 Service class TreadmillMeasures

Treadmill measures like speed, grade, time, distance, mets, heart rate, ...

Namespace: XDOP.schemas\_coscom\_org.service

Assembly: TreadmillMeasures.dll

### Syntax

```
public class TreadmillMeasures : IServiceObject
```

### Constructors

|  |  |
|--|--|
| <code>TreadmillMeasures(XDOPLib.IService service)</code> | Initializes a new instance of the class. |
|--|--|

### Common Service Properties

| Property name | Value  |
|---------------|--|
| ServiceId     | urn:schemas-coscom-org:serviceId:Measures          |
| MajorVersion  | 1  |
| MinorVersion  | 1  |
| ServiceType   | urn:schemas-coscom-org:service:TreadmillMeasures:1 |

### Service Properties and Events

| Type          | Name      | Flags | Description                              |
|---------------|-----------|-------|--|
| System.Single | Distance  | E G   | Actual distance value [m].               |
| System.Single | Energy    | E G   | Actual energy value [kJ].                |
| System.Single | Grade     | E G   | Actual grade value [%].                  |
| System.Byte   | HeartRate | E G   | Actual heart rate [1/min].               |
| System.Single | MET       | E G   | Actual MET (metabolic equivalent) value. |
| System.Single | Power     | E G   | Actual power value [W].                  |
| System.Single | Speed     | E G   | Actual speed value [m/s].                |
| System.UInt32 | Time      | E G   | Actual run time value [s].               |

### Service Methods

| Name                          | Description  |
|-------------------------------|--|
| <a href="#">ResetDistance</a> | Resets the distance value.                                   |
| <a href="#">ResetEnergy</a>   | Resets the energy value.                                     |
| <a href="#">ResetTime</a>     | Resets the run time.   |
| <a href="#">SetDistance</a>   | Sets the distance value.                                     |
| <a href="#">SetEventMask</a>  | Sets the event mask to control the eventing of each measure. |
| <a href="#">SetTime</a>       | Sets the time values.  |

### Service error codes

| Error code | Error description      |
|------------|------------------------|
| 830        | Distance value too big |

### Details

On the following pages you will find a detailed description for every action and variable of this service.

## 1.8.1 Actions

### 1.8.1.1 Action ResetDistance

Resets the distance value.

#### Syntax

```
public void ResetDistance()
```

### 1.8.1.2 Action ResetEnergy

Resets the energy value.

#### Syntax

```
public void ResetEnergy()
```

### 1.8.1.3 Action ResetTime

Resets the run time.

#### Syntax

```
public void ResetTime()
```

### 1.8.1.4 Action SetDistance

Sets the distance value.

#### Syntax

```
public void SetDistance(  
    System.Single distance  
)
```

#### Parameters

*distance*

The resolution of the distance value is 0.1 m and the maximum value depends on the treadmill type.

e.g. type 1.3 (167400 increments / m): 0 - 256569 m

#### Action Errors

|     |                        |
|-----|------------------------|
| 830 | Distance value too big |
|-----|------------------------|

### 1.8.1.5 Action SetEventMask

Sets the event mask to control the eventing of each measure.

## Syntax

```
public void SetEventMask(  
    System.Byte[] eventMask  
)
```

### Parameters

*eventMask*

Event Mask: time 0x01, distance 0x02, speed: 0x04, grade 0x08, heart rate 0x10, power 0x20, energy 0x40, MET 0x80

## 1.8.1.6 Action SetTime

Sets the time values.

## Syntax

```
public void SetTime(  
    System.Byte hrs  
    System.Byte min  
    System.Byte sec  
)
```

### Parameters

*hrs*

Hours value from 0 to 255

*min*

Minutes value:

hrs = 0: 0 - 59 minutes or 0 - 99 minutes

hrs > 0: 0 - 59 minutes

*sec*

Seconds value from 0 to 59

## 1.9 Service class MCU5FlashUpdate

This service can be used to update the firmware in the flash of the microcontroller. The update is separated in two phases. First a data buffer (data flash on the electronic board) is written with the firmware and constant data (WriteBuffer). In the second phase the microcontroller reads from the data flash part and writes to its on-chip flash memory (StartUpdate) without any intervention from the client. After a successful update of the on-chip flash, the microcontroller resets automatically.

Namespace: XDOP.schemas\_coscom\_org.service  
Assembly: MCU5FlashUpdate.dll

### Syntax

```
public class MCU5FlashUpdate : IServiceObject
```

### Constructors

|   |  |
|---|--|
| MCU5FlashUpdate(XDOPLib.IService service) | Initializes a new instance of the class. |
|---|--|

### Common Service Properties

| Property name | Value  |
|---------------|--|
| ServiceId     | urn:schemas-coscom-org:serviceId:MCU5FlashUpdate |
| MajorVersion  | 1  |
| MinorVersion  | 1  |
| ServiceType   | urn:schemas-coscom-org:service:MCU5FlashUpdate:1 |

### Service Properties and Events

| Type          | Name            | Flags | Description |
|---------------|-----------------|-------|-------------|
| System.UInt16 | UpdatePageCount | Edown |             |

### Service Methods

| Name                               | Description   |
|------------------------------------|---|
| <a href="#">GetDefaultType</a>     | Gets the default type and sub type, which has been set by the factory. This values are important, if the configuration parameters (options) of the treadmill are lost (e.g. the battery of the Real-Time-Clock has been replaced or the signature of the configuration memory is false).  |
| <a href="#">GetFirmwareVersion</a> | Gets the firmware version.  |
| <a href="#">ReadBuffer</a>         | Reads 128 bytes of the data flash buffer. The data flash buffer has a size of 256 bytes. An even page index reads the lower 128 bytes of data flash buffer and an odd page index reads the upper 128 bytes of the data flash buffer. Any call to ReadBuffer loads the data flash buffer from the data flash first. The data flash has a size of 256 kByte (= 2048 * 128). |
| <a href="#">ShowUpdateScreen</a>   | Shows flash update strings on the displays. This action should be used before the action "StartUpdate" is called.   |

|                             |  |
|-----------------------------|--|
| <a href="#">SoftReset</a>   | After the response of this action is sent, the microcontroller jumps to the reset address (soft reset).  |
| <a href="#">StartUpdate</a> | Starts the flash update of the microcontroller built-in flash.   |
| <a href="#">WriteBuffer</a> | Writes 128 bytes to the data flash buffer. The data flash buffer has a size of 256 bytes. If the second 128 Bytes (odd page index) of the buffer are written, then the whole buffer is written to the data flash. The data flash has a size of 256 kByte (= 2048 * 128). |

## Service error codes

No error codes defined

## Details

On the following pages you will find a detailed description for every action and variable of this service.

## 1.9.1 Actions

### 1.9.1.1 Action GetDefaultType

Gets the default type and sub type, which has been set by the factory. This values are important, if the configuration parameters (options) of the treadmill are lost (e.g. the battery of the Real-Time-Clock has been replaced or the signature of the configuration memory is false).

#### Syntax

```
public void GetDefaultType(  
    ref System.Byte type  
    ref System.Byte subType  
)
```

#### Parameters

*type*

Major type (major treadmill type)

*subType*

Sub type (treadmill type variant)

### 1.9.1.2 Action GetFirmwareVersion

Gets the firmware version.

#### Syntax

```
public System.String GetFirmwareVersion(  
    ref System.String version  
)
```

#### Parameters

*version*

String with 4 digits (e.g. 4011)

### 1.9.1.3 Action ReadBuffer

Reads 128 bytes of the data flash buffer. The data flash buffer has a size of 256 bytes. An even page index reads the lower 128 bytes of data flash buffer and an odd page index reads the upper 128 bytes of the data flash buffer. Any call to ReadBuffer loads the data flash buffer from the data flash first. The data flash has a size of 256 kByte (= 2048 \* 128).

#### Syntax

```
public void ReadBuffer(  
    System.UInt16 page128Idx  
    System.Byte key1  
    System.Byte key2  
    ref System.Byte[] packet  
)
```

#### Parameters

*page128Idx*

Index of the 128 byte long page:  
 $0 \leq \text{page128Idx} < 2048$

*key1*

Access key 1

*key2*

Access key 2

*packet*

Flash memory packet:  
128 flash buffer bytes + 1 checksum byte = 129 bytes  
129 bytes base64 coded = 172 characters  
 $(129 / 3) * 4 = 172$

#### 1.9.1.4 Action ShowUpdateScreen

Shows flash update strings on the displays. This action should be used before the action "StartUpdate" is called.

##### Syntax

```
public void ShowUpdateScreen()
```

#### 1.9.1.5 Action SoftReset

After the response of this action is sent, the microcontroller jumps to the reset address (soft reset).

##### Syntax

```
public void SoftReset()
```

#### 1.9.1.6 Action StartUpdate

Starts the flash update of the microcontroller built-in flash.

##### Syntax

```
public void StartUpdate(  
    System.Byte key1  
    System.Byte key2  
)
```

Parameters

*key1*

Access key 1



*key2*  
Access key 2

#### 1.9.1.7 Action WriteBuffer

Writes 128 bytes to the data flash buffer. The data flash buffer has a size of 256 bytes. If the second 128 Bytes (odd page index) of the buffer are written, then the whole buffer is written to the data flash. The data flash has a size of 256 kByte (= 2048 \* 128).

#### Syntax

```
public void WriteBuffer(  
    System.UInt16 page128Idx  
    System.Byte[] packet  
    System.Byte key1  
    System.Byte key2  
    ref System.Byte error  
)
```

#### Parameters

##### *page128Idx*

Index of the 128 byte long page:  
 $0 \leq \text{page128Idx} < 2048$

##### *packet*

Flash update packet:  
128 flash buffer bytes + 1 checksum byte = 129 bytes  
129 bytes base64 coded = 172 characters  
 $(129 / 3) * 4 = 172$

##### *key1*

Access key 1

##### *key2*

Access key 2

##### *error*

Packet checksum error:  
0 = okay  
1 = wrong

## 1.10 Service class DataLogger

This service provides a configurable data logger. The output of the logger is evented.

Namespace: XDOP.schemas\_coscom\_org.service

Assembly: DataLogger.dll

### Syntax

```
public class DataLogger : IServiceObject
```

### Constructors

|   |  |
|---|--|
| <code>DataLogger(XDOPLib.IService service)</code> | Initializes a new instance of the class. |
|---|--|

### Common Service Properties

| Property name | Value                                       |
|---------------|---|
| ServiceId     | urn:schemas-coscom-org:serviceId:Logger     |
| MajorVersion  | 1   |
| MinorVersion  | 1   |
| ServiceType   | urn:schemas-coscom-org:service:DataLogger:1 |

### Service Properties and Events

| Type          | Name   | Flags | Description                                     |
|---------------|--------|-------|---|
| System.String | Record | E G   | Record string with all configured column values |

### Service Methods

| Name                                 | Description                          |
|--------------------------------------|--------------------------------------|
| <a href="#">SetColumnDelimiter</a>   | Sets the column delimiter character. |
| <a href="#">SetRecordColumns</a>     | Sets the columns for the log record. |
| <a href="#">SetRecordEventPeriod</a> | Sets the event period.               |
| <a href="#">StartLogging</a>         | Starts the logging.                  |
| <a href="#">StopLogging</a>          | Stops the logging.                   |

### Service error codes

No error codes defined

### Details

On the following pages you will find a detailed description for every action and variable of this service.

### 1.10.1 Actions

#### 1.10.1.1 Action SetColumnDelimiter

Sets the column delimiter character.

##### Syntax

```
public void SetColumnDelimiter(  
    System.Char delimiter  
)
```

##### Parameters

*delimiter*

Delimiter character for delimiting the column values in the log record

##### Remarks

Comment for Action SetColumnDelimiter

#### 1.10.1.2 Action SetRecordColumns

Sets the columns for the log record.

##### Syntax

```
public void SetRecordColumns(  
    System.String columnList  
)
```

##### Parameters

*columnList*

String with semi-colon separated indexes of columns, which should be logged:

- 1 = Time
- 2 = Distance
- 3 = Speed
- 4 = Grade
- 5 = Heart Rate
- 6 = Power
- 7 = Energy
- 8 = MET

Example: 1;3;4 for time, speed and grade

#### 1.10.1.3 Action SetRecordEventPeriod

Sets the event period.

##### Syntax

```
public void SetRecordEventPeriod(  
    System.Byte period  
)
```

Parameters

*period*

Event period (time between two log record outputs) [s]

#### 1.10.1.4 Action StartLogging

Starts the logging.

[Syntax](#)

```
public void StartLogging()
```

#### 1.10.1.5 Action StopLogging

Stops the logging.

[Syntax](#)

```
public void StopLogging()
```